


## H - SERIES PHOTODIODE / RECEIVER



**Application Note** This unit is a high performance photodiode/receiver operated at ambient temperature with a dual gain FET input transimpedance amplifier. The output voltage is proportional to the input signal current:  $V_{out} = I_{signal} \cdot R_f$ . The PD/AMP is a DC coupled dual gain system. Care should be taken in shielding the unit from stray light during operation to prevent saturation of the amplifier (and potential failure).


<b>SPECIFICATIONS</b>	
<b>Detector Type</b>	<b>1.0 mm dia extended InGaAs Photodiode</b>
<b>Operating Temperature- °C</b>	<b>22</b>
<b>Operating Wavelength- μm</b>	<b>1.2 – 2.1</b>
<b>Responsivity- V/W @ pk</b>	<b><math>1.0 \times 10^7 / 10^6</math></b>
<b>Noise- V/Hz<sup>1/2</sup></b>	<b><math>2.5 \times 10^{-6} / 10^{-7}</math></b>
<b>NEP- W/Hz<sup>1/2</sup> @ pk</b>	<b><math>&lt; 2.5 \times 10^{-13}</math></b>
<b>Bandwidth (-3dB)- Hz</b>	<b>DC – 2kHz</b>
<b>Power Requirements</b>	<b>+/- 9 VDC to +/- 15 VDC</b>
<b>Connections</b>	<b>BNC signal output. Shielded power cable terminated with a DB-9 connector directly couples the unit with the PS -1 Low Noise Power Supply.</b>

**RoHS Compliant**



## H - SERIES PHOTODIODE / RECEIVER

### SPECTRAL RESPONSE - IGA1.9-series





---

### H - SERIES PHOTODIODE / RECEIVER

---

#### **OPERATING THE H-SERIES PHOTODIODE/AMPLIFIER**

**POWER SUPPLY:** A bipolar power supply is required, +,- 6VDC to +,-15VDC, 20mA. This means a +V, central/common ground and a -V connection - 3 wires total, to pins 6, 7, & 8 on the D-sub connector. The power supply pins should be bypassed physically close to the amplifier module. Double check wiring prior to turning on power. Improper /reverse wiring will damage the unit.

**GAIN SELECT:** The unit is supplied with a switch which provides a 10:1 HI/LO gain function. "UP" position is HI; "DOWN" position is LO gain. Consult the individual data sheet for specific values. The adjustable gain units have a single-turn control potentiometer which adds another x1 to x10 variable gain following the first stage. Clockwise rotation of the pot increases the gain.

**AMBIENT LIGHT:** Because of the high gains involved, the unit must be shielded from ambient background light during operation. Measurement errors and/or saturation can result from improper shielding.

**OUTPUT CONNECTION:** The signal output is thru a BNC connector (or BNC terminated cable in the case of the 2-color and adjustable gain units) located on the back of the module.